

Small Bowel Capsule Endoscopy



Your capsule endoscopy is scheduled for:

Additional instructions:

To ensure the best results for your capsule endoscopy of the small bowel, follow your doctor's instructions carefully.

Copyright ©2001-2007 Given Imaging Ltd. GIVEN, GIVEN & Design, PILLCAM, PILLCAM & Logo, PILLCAM IMAGING CAPSULE & Design, AGILE, RAPID, RAPID ACCESS, ORDERWIN, ORDER WHEN I NEED, FINGERS HOLDING A CAPSULE & Logo, FINGERS HOLDING PILLCAM CAPSULE & Logo, ICCE, ICCE Logos, and International Conference on Capsule Endoscopy are Trademarks and/or Registered Trademarks of Given Imaging Ltd. its subsidiaries and/or affiliates in the United States and/or other countries. All other company or product names are the trademarks or registered trademarks of their respective holders. All rights not expressly granted are reserved.

GMB-0237-01



Small Bowel Capsule Endoscopy




For more information please visit www.givenimaging.com.

How Does the Capsule Get Eliminated, and Will I Feel It Come Out?

The capsule is disposable and passes naturally with your bowel movement. You should not feel any pain or discomfort.

What are the Possible Complications of Capsule Endoscopy?

Although complications may occur, they are rare when doctors who are specially trained and experienced in this procedure perform the test. A potential risk could be retention of the capsule. It is important for you to recognize early signs of possible complications. If you have a fever after the test, trouble swallowing, or increasing chest or abdominal pain, contact your doctor immediately.





Your doctor may determine that a capsule endoscopy of the small bowel is necessary for further evaluation of your condition. Capsule endoscopy with the PillCam® SB video capsule will provide your doctor with pictures of your small intestine. Read this brochure to learn about the procedure and get answers to commonly asked questions.

Feel free to ask your doctor or nurse any additional questions you may have.

What is PillCam Capsule Endoscopy of the Small Bowel?

Capsule endoscopy with the PillCam SB video capsule enables your doctor to examine




your entire small intestine. Your doctor will have you ingest a video capsule that has its own camera and light source. You can move freely during the exam, which lasts about eight hours. While the video capsule travels through your body, it sends images to a data recorder you will wear close to your waist. Afterwards your doctor will view the images on a video monitor.

How Should I Prepare for the Procedure?

You will receive preparation instructions before the examination. An empty stomach allows optimal viewing conditions, so you should start a liquid diet after lunch the day prior to the examination and have nothing to eat or drink, including water, for approximately 10 hours before a small bowel examination. Your doctor will tell you when to start fasting. Tell your doctor in advance about any medications you take as you might need to adjust your usual dose for the examination. Inform your doctor if you have a pacemaker or other electromedical devices. Report any previous abdominal surgery, swallowing problem, or previous history of obstructions in the gastrointestinal tract.

What Can I Expect During PillCam Capsule Endoscopy?

The doctor or nurse will prepare you for the examination by applying a sensor array to your abdomen with adhesive sleeves. The video capsule is swallowed with water and passes naturally through your digestive tract while transmitting video images to a data recorder that you will wear close to your waist. You will be able to eat four hours after the capsule ingestion unless your doctor instructs you otherwise. The examination is over approximately eight hours after the capsule was ingested.



What Happens After PillCam Capsule Endoscopy?

At the end of the procedure, the data recorder and sensor array will be removed and the images acquired during your exam will then be downloaded to a workstation for your doctor to review. After ingesting the capsule and until it is excreted, you should not have a Magnetic Resonance Imaging (MRI) examination or be in the vicinity of an MRI machine.



How Will I Know the Results of the Procedure?

After you return the equipment, your doctor will process the information from the data recorder and will view a color video of the pictures taken by the capsule. After the doctor has looked at this video, you will be contacted with the results.

IMPORTANT REMINDER: This brochure provides general information only and should not be used as a definitive basis for diagnosis or treatment in any particular case. It is very important that you consult your doctor about your specific condition.